



Hillrom offers one of the largest ranges of slings and lifting accessories. Our slings can be used to lift every unique person in all kinds of situations. The right choice of sizes is essential for safety, function and comfort. We offer a wide range of sizes for all our slings.

Designation	Article No.	Size	Weight Class	No cleaning needed/ possible	Sling could be left underneath the patient	seated position, so patient hip bend is required	extra head support for people with no/less head control	This product is already in-house, please tick	
SLINGS FOR TRANSFER OF PATIENTS IN UPRIGHT SITTING POSITION									
	Modell 200 Highback NET POLYESTER	35200115	M	40-80 kg max. 200 kg				<input type="checkbox"/>	
		35200116	L	70-120 kg max. 200 kg			✓	✓	<input type="checkbox"/>
		35200117	XL	120-250 kg max. 500kg					<input type="checkbox"/>
	Modell 200 Highback NET POLYESTER ■ suitable for transfer in most environments ■ Ideal for people with poor torso stability ■ Allows semi-reclining seated posture for optimum comfort	35200315	M	40-80 kg max. 200 kg				<input type="checkbox"/>	
		35200316	L	70-120 kg max. 200 kg		✓ if back support is removed	✓	✓	<input type="checkbox"/>
		35200317	XL	120-250 kg max. 500kg					<input type="checkbox"/>
	Model 25 Highback COTTON POLYESTER ■ disposable product, yet strong	3525915-2	M	40-80 kg max. 200 kg				<input type="checkbox"/>	
		3525916-2	L	70-120 kg max. 200 kg	✓		✓	✓	<input type="checkbox"/>
		3525917-2	XL	120-250 kg max. 500kg					<input type="checkbox"/>
	Model 000 POLYESTER	35000115	M	50-90 kg max. 200 kg				<input type="checkbox"/>	
		35000116	L	90-110 kg max. 200 kg			✓		<input type="checkbox"/>
		35000117	XL	<110 kg max. 500kg					<input type="checkbox"/>
	Hygiene Model 41 POLYESTER	3541135	M	60-80 kg max. 200 kg				<input type="checkbox"/>	
		3541136	L	80-100 kg max. 200 kg					<input type="checkbox"/>
	Model 000 NET POLYESTER	35000315	M	60-80 kg max. 200 kg				<input type="checkbox"/>	
		35000316	L	80-100 kg max. 200 kg				<input type="checkbox"/>	
		35000317	XL	100-150 kg max. 300kg				<input type="checkbox"/>	
	Sliding sheet short, no handles	3710015	M			dimensions: 70 x 90 cm		<input type="checkbox"/>	
SLINGS FOR HORIZONTAL TRANSFER AND MOBILISATION									
	Solo Reposheet COTTON/POLYESTER ■ disposable product, yet strong	36871004	XL	max. 300 kg	✓	✓	for repositioning of patients in bed, horizontal transfer with Ultra- and Flexo stretcher is also possible	<input type="checkbox"/>	
	Octostretch	3683106	L	max. 200 kg			use only together with Octostretch stretcher	<input type="checkbox"/>	
	Liftsheet	3684106	XL	max. 400 kg			use only together with Ultra- or Flexo stretcher	<input type="checkbox"/>	
	Solo Octosheet COTTON/POLYESTER ■ disposable product, yet strong	3683906-4	L	max. 200 kg	✓		use only together with Octostretch stretcher	<input type="checkbox"/>	