

Affinity® 4 Birthing Bed



Delivering Safety, Comfort and Convenience

The labor and delivery experience can be unpredictable for everyone involved. To help manage this, Hill-Rom offers a solution that can adapt and deliver for mom, baby and you.

The Affinity® 4 Birthing Bed is designed to support everyone throughout the birthing process.

Delivering Safety

- Infinite Patient Positioning
- Patient Transfers
- Patient Transports

Delivering Comfort

- Automatic Pelvic Tilt
- Wide Surface
- Lumbar Support

Delivering Convenience

- Battery Backup
- Clinical Access
- Infinite Positioning Calf Supports

Every mom is different, every labor is different. The Affinity 4 Birthing Bed is designed to help you adapt by delivering **safety, comfort and convenience**.

Delivering Safety, Comfort and Convenience



PATIENT TRANSPORT
Ergonomic push handles and light-weight birthing bed frame.

LUMBAR & SEAT SUPPORT

Customize patient comfort with the push of a button.

MAXIMUM PATIENT WEIGHT

500 lb bed patient weight limit.
400 lb foot section patient weight limit.



BATTERY BACKUP
Enables function of the bed during transport or loss of power.

GAS ASSIST TRENDELENBURG

Achieve 8° of Trendelenburg on demand.



15° PELVIC TILT
Tilted seat designed to cradle mom to increase her feeling of security.

PATIENT TRANSFER

Zero-transfer gap for a seamless transfer of your patient.

"DUAL-LOCKING" CASTERS

Two points of contact prevents casters from rolling or swivelling.



CLINICAL ACCESS
Available in both traditional U-Cut or our exclusive V-Cut mattress to improve perineal access.

STOW & GO® FOOT SECTION

Exclusive, ergonomic design; foot section slides and stores under the bed.



LIGHT-WEIGHT FOOT SECTION
Weighs 16 lbs, includes ergonomic handles and angular stand to keep mattress off the floor.

WIDE SLEEP DECK

35" surface width gives mom room to care for baby in bed.



INFINITE POSITIONING CALF SUPPORTS
Ergonomically designed, provides comfort and infinite positioning.

DELIVERING SAFETY

The Affinity® 4 Birthing Bed is designed to keep mom, baby and you safe.

DELIVERING COMFORT

Keeping mom comfortable throughout the birthing experience should be paramount.

DELIVERING CONVENIENCE

The Affinity® 4 Birthing Bed is designed to make every movement fast and easy.

Affinity® 4 Birthing Bed

Delivering Safety, Comfort and Convenience

Technical Specifications

| | |
|---|--|
| Product Number | P3700 |
| Full Bed Length | 90" (229 cm) P3700A-P3700C 92" (234 cm) P3700D-P3700E |
| Bed Width | |
| Siderails Up | 39" (99 cm) |
| Siderails Down | 36" (91 cm) |
| Bed Height | |
| Low (with mattress) | 24" (61 cm) |
| High (with mattress) | 39" (99 cm) |
| Maximum height of seat section in Trend-like position | 40" (102 cm) |
| Mattress | |
| Length | 78" (198 cm) |
| Width | 35" (89 cm) |
| Thickness (head section) | 5" (12.7 cm) |
| Thickness (foot section) | 3.5" (8.9 cm) |
| Critical Angles | |
| Maximum head elevation | 63° |
| Maximum seat elevation | 15° |
| Maximum Trend-like positioning | 8° |
| Maximum Lift Capacity | |
| Bed | 500 lbs (227 kg) |
| Foot section | 400 lbs (181 kg) |
| Caster size | 6" (15 cm) |
| Bed weight | 480 lbs (218 kg) |

Standard Features

V-Cut or Straight Cut Surfaces
Removable Foot Section
Battery Backup
Powered Foot Section
Four wheel dual-locking casters
OneStep® Siderail Release
Plastic Blow-Molded Head Panel
Automatic Nightlight
Lock-out controls
Manual Trend-like positioning
Instant CPR
Instant Labor Grips
Automatic Seat Tilt (15°)
Siderail Controls
IV Pole
Plastic Blow Molded Head Panel
Hinged Foot Section Mattress
Central Brake and Steer
Fluid Basin

Optional Features

Prevention Surface
Air Package
Nurse Call, Lighting, Entertainment & UTV
Contemporary Headboard
Freedom Hill™ Headboard
Liberty Hill™ Headboard
Calf Supports
Push Handles with IV Hooks

Accessories

Universal Instrument Tray
Labor Bar
Drainage Bag Hook Kit
Permanent IV Pole
Line Managers
Foot Rest Bar (for lift off beds)
Pendant Holder
Pendant Kit
Full Leg Supports
Full Leg Supports, Telescoping
Drapes - Medline/Proxima
Anesthesia Screen
Oxygen Tank Holder

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
207598 rev 1 12-MAR-2018 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Enhancing outcomes for patients and their caregivers:

Hill-Rom